

Software Releases and Updates

SLIM			
	Version	Release Date	Description
	3.4.0	09 Nov 2006	Spreadsheet Data Entry is supported and other significant improvements.
	2.12.0	09 Nov 2006	Spreadsheet Data Entry is supported and other significant improvements
	SP03080612A	25 Aug 2008	Formulation names are no longer required to be unique.
	SP03080129A	01 Feb 2008	Display of very large reports.
	SP03071105A	21 Dec 2007	Performance Improvements.
	SP03070619A	15 Aug 2007	LDAP Connectivity added for Windows Server 2003.
	SP03070219A	16 Apr 2007	Text specifications were added to reports.
	SP03061130A	03 Jan 2007	Corrected defects associated with version 3.4.0
SLIM Report Generator			
	Version	Release Date	Description
	2.2.0 (SQL)	06 Oct 2006	Many configurable features were added.
	2.2.0 (Oracle)	01 Dec 2006	Corrected minor defects associated with version 2.1.0
	SP03A071005A	21 Dec 2007	Added packaging description and fixed some international settings.
SLIMStat			
	Version	Release Date	Description
	4.3.0	15 Jan 2008	Vista Compatible and other new features.
	4.2.0	27 Jun 2005	Model 3 and Model 4 date pooling is supported and other significant improvements.
	SP04060308A	06 Apr 2006	Prevent File Save As from being displayed when File Save is selected.

SLIMStat (continued)			
	Version	Release Date	Description
	SP04051214A	18 Jan 2006	The Model 2 pooling confidence interval calculation has been changed so that under certain circumstances the expiration date will yield more conservative values.
SLIM (Oracle)			
	Version	Release Date	Description
	3.5.1	10 Nov 2008	Minor fixes to Server version 3.5.0.
	2.13.1	10 Nov 2008	Minor fixes to Client version 2.13.0.
	3.5.0	01 Oct 2008	Additional Formulation functionality is supported and other significant improvements.
	2.13.0	01 Oct 2008	Additional Formulation functionality is supported and other significant improvements.
	3.4.0	09 Mar 2006	Spreadsheet Data Entry is supported and other significant improvements.
	2.11.0	09 Mar 2006	Spreadsheet Data Entry is supported and other significant improvements.
	SP07081202A	17 Dec 2008	Correct procedure pertaining to email notification.
	SP07080718A	23 Oct 2008	Fixes to the Result Entry Spreadsheet.
	SP07060614A	07 Jul 2006	Fixes to the stability study framework report.
	SP07060309A	17 May 2007	Support Oracle Client Version 10.2.0.1.
	SP07060327A	11 Apr 2006	Fixes to the data entry spreadsheet.

PC/Chrom			
	Version	Release Date	Description
	4.5.0	01 Oct 2008	TCP/IP (Ethernet) A/D Hardware is now supported. The default A/D input voltage range was increased from +/- 1250 mV to +/- 2000 mV. A data acquisition rate of 32 Hz is now supported.
	4.4.0	14 Jan 2008	The automatic peak integration algorithm has been optimized to prevent integration of negative areas of chromatographic peaks and the Baseline Noise algorithm has been optimized so that reasonable changes in the slope of the baseline do not adversely impact the signal-to-noise (S/N) calculation.
	4.3.1	17 May 2006	Corrected minor defects associated with version 4.3.0
ChromWave			
	Version	Release Date	Description
	3.1.0	24 Apr 2006	Added improvements.
	3.0.0	29 Nov 2005	Initial Release with 21 CFR Part compliance features.
IntelliFORM			
	Version	Release Date	Description
	1.0.0	29 Mar 2001	Initial Release
IntelliLIMS			
	Version	Release Date	Description
	1.2.0	08 May 2008	Subsamples can now be extracted from existing samples with a separate chain of custody.
	1.1.0	01 Oct 2007	Minor updates and improvements.

IntelliLIMS (continued)			
	Version	Release Date	Description
	SP09080317A	14 Jan 2009	Improved Sample Manager performance, as well as other minor updates and improvements.
	1.0.0	30 Jul 2007	Initial Release